


## Anti-Phospho-Tbk1/Nak-S172 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Uhd2

Item Number: CM0000706



### Introduction

High-quality rabbit monoclonal antibody specific for phosphorylated TBK1 (Ser172). Validated for WB and ELISA, cross-reacts with human and mouse. Ideal for innate immunity, autophagy, and inflammation research.

[Learn More](#)

| Parameter        | Value                                                            |
|------------------|------------------------------------------------------------------|
| Product Name     | Phospho-TBK1/NAK-S172 Rabbit mAb                                 |
| Remarks/Alias    | NAK; T2K; IIAE8; FTDALS4; Phospho-TBK1/NAK-S172                  |
| Species          | Human                                                            |
| Gene ID          | 29110 (Human)                                                    |
| Immunogen        | A synthetic phosphorylated peptide around S172 of human TBK1/NAK |
| Source           | Rabbit                                                           |
| Category         | Phosphorylated Antibodies                                        |
| Application      | WB, ELISA                                                        |
| Cross Reactivity | Human, Mouse                                                     |
| SWISS            | Q9UHD2                                                           |
| Protein Weight   | 84kDa                                                            |
| Shipping         | Ice bag                                                          |