


Anti-Phospho-Tbk1/Nak-S172 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uhd2

Item Number: CM0001369



Introduction

This rabbit polyclonal antibody detects phosphorylated TBK1/NAK at Ser172. Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Suitable for innate immunity, antiviral signaling, and autophagy studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TBK1/NAK-S172 Rabbit pAb
Remarks/Alias	NAK; T2K; IIAE8; FTDALS4; Phospho-TBK1/NAK-S172
Species	Human
Gene ID (Human)	29110
Gene ID	29110
Immunogen	A synthetic phosphorylated peptide around S172 of human TBK1/NAK (NP_037386.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHD2
Protein Weight	84 kDa
Shipping	Ice bag