

Anti-Phospho-Tdp-43/Tardbp (S409/S410) Rabbit Polyclonal Antibody For Wb, Db, Elisa - Q13148

Item Number: CM0001763



Introduction

Phospho-specific rabbit polyclonal antibody targeting S409/S410 of human TDP-43 (TARDBP). Validated for WB, DB, and ELISA. Suitable for studying neurodegenerative disease mechanisms and RNA biology. Detects a 45 kDa band.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TDP-43/TARDBP-S409/S410 Rabbit pAb
Remarks/Alias	ALS10; TDP-43; Phospho-TDP-43/TARDBP-S409/S410
Species	Human
GeneID (Human)	23435
GeneID	23435
Immunogen	A synthetic phosphorylated peptide around S409&S410 of human TDP-43/TARDB (NP_031401.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human
SWISS	Q13148
Protein Weight	45kDa
Shipping	Ice bag