


Anti-Phospho-Tifa-T9 Rabbit Monoclonal Antibody For Wb, Elisa - Q96Cg3

Item Number: CM0001178



Introduction

Phospho-TIFA-T9 Rabbit mAb detects phosphorylated TIFA at Thr-9 in human samples. Suitable for WB and ELISA, ideal for innate immunity and NF-κB signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TIFA-T9 Rabbit mAb
Remarks/Alias	T2BP; T6BP; TIFAA; Phospho-TIFA-T9
Species	Human
GeneID (Human)	92610
GeneID	92610
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96CG3
Protein Weight	21kDa
Shipping	Ice bag