


Anti-Phospho-Trka-Y490 Rabbit Polyclonal Antibody For Wb And Elisa - P04629

Item Number: CM0000840



Introduction

Phospho-specific rabbit polyclonal antibody against TrkA-Y490 (NTRK1) validated for Western blot and ELISA. Suitable for human, mouse, and rat samples in neurotrophic signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TrkA-Y490 Rabbit pAb
Remarks/Alias	MTC; TRK; TRK1; TRKA; Trk-A; p140-TrkA; Phospho-TrkA-Y490
Species	Human
GeneID (Human)	4914
GeneID	4914
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y490 of human TrkA (NP_001012331.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04629
Protein Weight	87kDa
Shipping	Ice bag