

Anti-Phospho-Trka-Y674/675 + Trkb-Y706/707 Rabbit Polyclonal Antibody For Wb, Elisa - P04629-2 / Q16620

Item Number: CM0001463

Introduction

Phospho-TrkA-Y674/675 + TrkB-Y706/707 Rabbit pAb detects phosphorylated human TrkA (Y674/675) and TrkB (Y706/707). Validated for WB and ELISA. Ideal for neurotrophin signaling and receptor activation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TrkA-Y674/675 + TrkB-Y706/707 Rabbit pAb
Remarks/Alias	NTRK1/NTRK2; Phospho-TrkA-Y674/675 + TrkB-Y706/707
Species	Human
Gene ID (Human)	4914/4915
Gene ID	4914/4915
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y674 & Y675 of human TrkA/NTRK1 (NP_001012331.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04629-2/Q16620
Protein Weight	35kDa/53-60kDa/77kDa/83-93kDa
Shipping	Ice bag