


Anti-Phospho-Tsc2 (S1387) Rabbit Polyclonal Antibody For Wb And Elisa - P49815

Item Number: CM0001821



Introduction

Phospho-specific rabbit polyclonal antibody against Tuberin (TSC2) phosphorylated at Ser1387. Validated for Western blot and ELISA with human, mouse, and rat reactivity. Ideal for mTOR pathway and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TSC2-S1387 Rabbit pAb
Remarks/Alias	LAM; TSC4; PPP1R160; Phospho-TSC2-S1387
Species	Human
GeneID (Human)	7249
GeneID	7249
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49815
Protein Weight	201kDa
Shipping	Ice bag