


Anti-Phospho-Tsc2-T1462 Rabbit Polyclonal Antibody For Wb - P49815

Item Number: CM0000696



Introduction

Phospho-specific rabbit polyclonal antibody detecting TSC2 (Tuberin) phosphorylation at T1462. Suitable for WB and ELISA in human and mouse. Ideal for studying mTORC1 regulation and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TSC2-T1462 Rabbit pAb
Remarks/Alias	LAM; TSC4; PPP1R160; Phospho-TSC2-T1462
Species	Human
Gene ID (Human)	7249
Gene ID	7249
Immunogen	A synthetic phosphorylated peptide around T1462 of human TSC2 (NP_000539.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49815
Protein Weight	201kDa
Shipping	Ice bag