

Anti-Phospho-Ubiquitin-Ser65 Rabbit Polyclonal Antibody For Wb, Elisa - P62979 / P62987 / P0Cg47 / P0Cg48

Item Number: CM0001301



Introduction

Anti-Phospho-Ubiquitin-Ser65 Rabbit Polyclonal Antibody validated for WB, ELISA. Detects phosphorylated Ser65 of human ubiquitin; ideal for studying ubiquitination, protein degradation, and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Ubiquitin-Ser65 Rabbit pAb
Remarks/Alias	UBC; S27A; CEP80; UBA80; HEL112; UBCEP1; UBCEP80; L40; CEP52; RPL40; HUBCEP52; HEL-S-50; HMG20; Phospho-Ubiquitin-Ser65
Species	Human
GeneID (Human)	6233/7311/7314/7316
GeneID	6233/7311/7314/7316
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S65 of human Ubiquitin (NP_002945.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P62979/P62987/P0CG47/P0CG48
Protein Weight	18kDa/15kDa/77kDa
Shipping	Ice bag