

## Anti-Phospho-B-Catenin-S45 Rabbit Polyclonal Antibody For Wb, Elisa - P35222

Item Number: CM0000912



### Introduction

Phospho-β-Catenin-S45 rabbit polyclonal antibody (P35222) for WB and ELISA. Detects endogenous phosphorylated β-catenin at Ser45 in human, mouse, and rat. Ideal for Wnt signaling and cancer research.

[Learn More](#)

| Parameter        | Value                                                                           |
|------------------|---------------------------------------------------------------------------------|
| Product Name     | Phospho-β-Catenin-S45 Rabbit pAb                                                |
| Remarks/Alias    | EVR7; CTNNB; MRD19; NEDSDV; armadillo; Phospho-β-Catenin-S45                    |
| Species          | Human                                                                           |
| GeneID (Human)   | 1499                                                                            |
| GeneID           | 1499                                                                            |
| Immunogen        | A synthetic phosphorylated peptide around S45 of human β-catenin (NP_001091679) |
| Source           | Rabbit                                                                          |
| Category         | Phosphorylated Antibodies                                                       |
| Application      | WB, ELISA                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                               |
| SWISS            | P35222                                                                          |
| Protein Weight   | 85kDa                                                                           |
| Shipping         | Ice bag                                                                         |