

Anti-Phospho-B-Catenin (S552) Polyclonal Antibody For Wb, Elisa - P35222

Item Number: CM0000871



Introduction

Rabbit polyclonal antibody targeting phospho-β-Catenin (S552). Validated for WB, ELISA with human, mouse, rat cross-reactivity. Ideal for Wnt signaling and cell adhesion research. - P35222

[Learn More](#)

Parameter	Value
Product Name	Phospho-β-Catenin-S552 Rabbit pAb
Remarks/Alias	EV7; CTNNB; MRD19; NEDSDV; armadillo; Phospho-β-Catenin-S552
Species	Human
GeneID (Human)	1499
GeneID	1499
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S552 of human beta Catenin (NP_001091679)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35222
Protein Weight	85 kDa
Shipping	Ice bag