

Anti-Phospho-Vangl2 (S79/S82/S84) Rabbit Monoclonal Antibody For Wb And Elisa - Q9UIk5

Item Number: CM0001294



Introduction

Rabbit monoclonal antibody specific to phosphorylated VANGL2 at serines 79, 82, and 84. Validated for Western blot and ELISA, reacts with human and rat samples. Suitable for studying planar cell polarity and developmental processes.

[Learn More](#)

Parameter	Value
Product Name	Phospho-VANGL2-S79/S82/S84 Rabbit mAb
Remarks/Alias	LPP1; LTAP; STB1; STBM; STBM1; Phospho-VANGL2-S79/S82/S84
Species	Human
GeneID (Human)	57216
GeneID	57216
Immunogen	Synthetic Peptide: A phospho specific peptide corresponding to residues surrounding S79/S82/S84 of human VANGL2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9ULK5
Protein Weight	60kDa
Shipping	Ice bag