

Anti-Phospho-Vasp (S157) Rabbit Monoclonal Antibody For Wb, Elisa - P50552

Item Number: CM0001121



Introduction

Anti-Phospho-VASP (S157) Rabbit Monoclonal Antibody validated for WB and ELISA in human samples. Targets phosphorylated VASP at S157, an actin-associated protein pivotal in cytoskeleton remodeling, platelet aggregation, and cell motility. Suitable for vascular biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-VASP-S157 Rabbit mAb
Remarks/Alias	VASP; Phospho-VASP-S157
Species	Human
Gene ID (Human)	7408
Gene ID	7408
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S157 of human VASP (NP_003361.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P50552
Protein Weight	40kDa
Shipping	Ice bag