

Anti-Phospho-Vasp-S157 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P50552

Item Number: CM0000791



Introduction

Phospho-specific rabbit polyclonal antibody targeting VASP (phospho-S157). Validated for WB, IHC-P, IP, ELISA with human, mouse, rat reactivity. VASP regulates actin dynamics, platelet aggregation, and cell migration. Ideal for cytoskeleton and thrombosis research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-VASP-S157 Rabbit pAb
Remarks/Alias	VASP; Phospho-VASP-S157
Species	Human
GeneID (human)	7408
GeneID	7408
Immunogen	Synthetic phosphorylated peptide around S157 of human VASP (NP_003361.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50552
Protein Weight	40kDa
Shipping	Ice bag