

Anti-Phospho-Vasp-S239 Polyclonal Antibody For Wb, Ihc-P, Elisa - P50552

Item Number: CM0000769



Introduction

Rabbit polyclonal antibody detects Phospho-VASP-S239 (Vasodilator-stimulated phosphoprotein). Validated for Western Blot, IHC-P, ELISA. Reacts with Human, Mouse, Rat. Ideal for cytoskeleton and cell migration research.

[Learn More](#)

Specification	Details
Product Name	Phospho-VASP-S239 Rabbit pAb
Remarks/Alias	VASP; Phospho-VASP-S239
Species	Human
GeneID (human)	7408
GeneID	7408
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S239 of human VASP (NP_003361.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50552
Protein Weight	40kDa
Shipping	Ice bag