

Anti-Phospho-Vegf Receptor 2-Y1175 Rabbit Polyclonal Antibody For Wb, Elisa - P35968

Item Number: CM0000760



Introduction

Rabbit polyclonal antibody targeting phosphorylated Y1175 of VEGFR2/KDR. Validated for Western blot and ELISA, cross-reacts with human, mouse, and rat. Ideal for angiogenesis and cancer signaling research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-VEGF Receptor 2-Y1175 Rabbit pAb
Remarks/Alias	FLK1, CD309, VEGFR, VEGFR2, Phospho-VEGF Receptor 2-Y1175
Species	Human
Gene ID	3791
Immunogen	A synthetic phosphorylated peptide around Y1175 of human VEGF Receptor 2 (NP_002244.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35968
Protein Weight	152kDa
Shipping	Ice bag