

Anti-Phospho-Yap1-S127 Rabbit Monoclonal Antibody For Wb, Elisa - P46937

Item Number: CM0001733



Introduction

High-quality rabbit monoclonal antibody specific to phosphorylated YAP1 at serine 127. Validated for Western blot and ELISA in human, mouse, and rat samples. Ideal for Hippo signaling pathway and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-YAP1-S127 Rabbit mAb
Remarks/Alias	YAP1; COB1; YAP; YAP2; YAP65; YKI; Yes associated protein 1; Phospho-YAP1-S127
Species	Human
GeneID (human)	10413
GeneID	10413
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S127 of human YAP1 (NP_001123617.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46937
Protein Weight	36kDa/48kDa/49kDa/50kDa/52kDa/53kDa/54kDa
Shipping	Ice bag