

## Anti-Phospho-Yap1-S127 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P46937

Item Number: CM0000970



### Introduction

Phospho-specific YAP1 (Ser127) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets the Hippo pathway transcriptional coactivator YAP1, a key regulator of organ size and tumor suppression. Suitable for cancer and signaling research.

[Learn More](#)

| Parameter        | Value                                                                       |
|------------------|-----------------------------------------------------------------------------|
| Product Name     | Phospho-YAP1-S127 Rabbit pAb                                                |
| Remarks/Alias    | YAP; YKI; COB1; YAP2; YAP-1; YAP65; Phospho-YAP1-S127                       |
| Species          | Human                                                                       |
| GeneID (Human)   | 10413                                                                       |
| GeneID           | 10413                                                                       |
| Immunogen        | Synthetic phosphorylated peptide around S127 of human YAP1 (NP_001123617.1) |
| Source           | Rabbit                                                                      |
| Category         | Phosphorylated Antibodies                                                   |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                           |
| SWISS            | P46937                                                                      |
| Protein Weight   | 54kDa                                                                       |
| Shipping         | Ice bag                                                                     |