


Anti-Phospho-Yap1-S128 Rabbit Polyclonal Antibody For Wb - P46937

Item Number: CM0001469



Introduction

Phospho-specific rabbit polyclonal antibody targeting YAP1 phosphorylated at Ser128, a key Hippo pathway effector. Reacts with human and mouse samples, validated for Western blot and ELISA applications. Ideal for cancer and organ size research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-YAP1-S128 Rabbit pAb
Remarks/Alias	YAP; YKI; COB1; YAP2; YAP-1; YAP65; Phospho-YAP1-S128
Species	Human
Gene ID (Human)	10413
Gene ID	10413
Immunogen	A synthetic phosphorylated peptide around S128 of human YAP1 (NP_001123617.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46937
Protein Weight	54kDa
Shipping	Ice bag