

Anti-Phospho-Ywhaz-S58 Rabbit Pab For Wb And If/Icc - P63104

Item Number: CM0001098



Introduction

Rabbit polyclonal antibody specifically detects YWHAZ phosphorylated at serine-58. Suitable for western blot and immunofluorescence in human, mouse, and rat samples. Ideal for signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-YWHAZ-S58 Rabbit pAb
Remarks/Alias	HEL4; YWHAD; KCIP-1; HEL-S-3; POPCHAS; HEL-S-93; 14-3-3-zeta; Phospho-YWHAZ-S58
Species	Human
Gene ID	7534
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human, Mouse, Rat
SWISS	P63104
Protein Weight	28kDa
Shipping	Ice bag