


Anti-Phospho-Zap70-Y319 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P43403

Item Number: CM0001150



Introduction

Phospho-ZAP70-Y319 rabbit polyclonal antibody (pAb) validated for Western blotting (WB), immunoprecipitation (IP), and ELISA in human samples. Specific for the phosphorylated Y319 residue of ZAP70, a key kinase in T-cell receptor signaling and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ZAP70-Y319 Rabbit pAb
Remarks/Alias	SRK; STD; TZK; STCD; IMD48; ADMIO2; ZAP-70; Phospho-ZAP70-Y319
Species	Human
GeneID (Human)	7535
Immunogen	A synthetic phosphorylated peptide around Y319 of human ZAP70 (NP_001070.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P43403
Protein Weight	70kDa
Shipping	Ice bag