


Anti-Phospho-Zap70-Y319+Syk-Y352 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P43403/P43405

Item Number: CM0001591



Introduction

Rabbit polyclonal antibody recognizing phospho-Zap70 (Y319) and phospho-Syk (Y352). Suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Useful for T-cell and B-cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Zap70-Y319+Syk-Y352 Rabbit pAb
Remarks/Alias	SRK; STD; TZK; STCD; IMD48; ADMIO2; ZAP-70
Species	Human
GeneID (Human)	7535
GeneID	7535/6850
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y319 of human ZAP70 (NP_001070.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P43403/P43405
Protein Weight	70kDa
Shipping	Ice bag