


Anti-Phospho-Zap70 (Y493) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P43403

Item Number: CM0001417



Introduction

Phospho-specific rabbit polyclonal antibody targeting ZAP70 at tyrosine 493, validated in WB, IHC-P, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Ideal for T-cell receptor signaling and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ZAP70-Y493 Rabbit pAb
Remarks/Alias	SRK; STD; TZK; STCD; IMD48; ADMIO2; ZAP-70; Phospho-ZAP70-Y493
Species	Human
Gene ID (Human)	7535
Gene ID	7535
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y493 of human ZAP70 (NP_001070.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43403
Protein Weight	70kDa
Shipping	Ice bag