


Anti-Phospholamban Monoclonal Antibody For Wb, If-P, Elisa - P26678

Item Number: CM0007076



Introduction

Anti-Phospholamban rabbit monoclonal antibody validated for WB, IF-P, ELISA with cross-reactivity to mouse and rat. Targets the sarcoplasmic reticulum calcium regulator involved in cardiac muscle contractility, essential for cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	Phospholamban Rabbit mAb
Remarks/Alias	PLN; CMD1P; CMH18; PLB; phospholamban; Phospholamban
Species	Human
Gene ID	5350
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-52 of human Phospholamban. (NP_002658.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P26678
Protein Weight	6 kDa
Shipping	Ice bag