


Anti-Phospholamban (PLN) Polyclonal Antibody For Wb, Elisa - P26678

Item Number: CM0021750



Introduction

Rabbit polyclonal antibody against human Phospholamban (PLN) for Western blot and ELISA, with cross-reactivity to mouse and rat. Targets a key regulator of cardiac calcium homeostasis and contractility.

[Learn More](#)

Parameter	Value
Product Name	PLN Rabbit pAb
Remarks/Alias	PLB; CMD1P; CMH18; PLN
Species	Human
GeneID (Human)	5350
GeneID	5350
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-52 of human PLN (NP_002658.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P26678
Protein Weight	6 kDa
Shipping	Ice bag