

Anti-Phospholipid Scramblase 1 (Plscr1) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15162

Item Number: CM0008775



Introduction

Phospholipid Scramblase 1 (PLSCR1) Rabbit mAb validated for WB, IHC-P, IF/ICC, ELISA with reactivity across human, mouse, and rat—ideal for apoptosis and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	Phospholipid Scramblase 1 (PLSCR1) Rabbit mAb
Remarks/Alias	MMTRA1B; Phospholipid Scramblase 1 (PLSCR1)
Species	Human
GeneID (Human)	5359
GeneID	5359
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Phospholipid Scramblase 1 (PLSCR1) (NP_066928.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15162
Protein Weight	35kDa
Shipping	Ice bag