


Anti-Phospholipid Scramblase 1 (Plscr1) Polyclonal Antibody For Wb And Elisa - O15162

Item Number: CM0031052



Introduction

Rabbit polyclonal antibody against human Phospholipid Scramblase 1 (PLSCR1). Validated for WB and ELISA; detects 35kDa protein. Ideal for apoptosis, viral defense, and lipid scrambling research.

[Learn More](#)

Product Name	Phospholipid Scramblase 1 (PLSCR1) Rabbit pAb
Remarks/Alias	MMTRA1B; Phospholipid Scramblase 1 (PLSCR1)
Species	Human
GeneID (human)	5359
GeneID	5359
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-318 of human Phospholipid Scramblase 1 (PLSCR1) (NP_066928.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15162
Protein Weight	35kDa
Shipping	Ice bag