

Anti-Phox2B Rabbit Monoclonal Antibody For Wb, Elisa - Q99453

Item Number: CM0012625



Introduction

Rabbit monoclonal antibody against human PHOX2B, validated for WB and ELISA. Cross-reacts with mouse. Targets a transcription factor crucial for noradrenergic neuron development and implicated in neuroblastoma and CCHS. Suitable for neuroscience and oncology research.

[Learn More](#)

Parameter	Specification
Product Name	PHOX2B Rabbit mAb
Remarks/Alias	CCHS; PMX2B; NBLST2; NBPhox; PHOX2B
Species	Human
Gene ID (Human)	8929
Gene ID	8929
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 215-314 of human PHOX2B (NP_003915.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99453
Protein Weight	32kDa
Shipping	Ice bag