


Anti-Phox2B Rabbit Polyclonal Antibody For Wb And Elisa - Q99453

Item Number: CM0027982



Introduction

Anti-PHOX2B Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human, mouse, rat PHOX2B, a transcription factor essential for autonomic neuron development. Ideal for neurobiology, neuroblastoma, and developmental disease research.

[Learn More](#)

Parameter	Value
Product Name	PHOX2B Rabbit pAb
Remarks/Alias	CCHS, PMX2B, NBLST2, NBPhox, PHOX2B
Species	Human
Gene ID (Human)	8929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human PHOX2B (NP_003915.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99453
Protein Weight	32kDa
Shipping	Ice bag