

Anti-Phpt1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nrx4

Item Number: CM0013469



Introduction

Rabbit polyclonal antibody targeting human PHPT1 (14 kDa phosphohistidine phosphatase). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for research on histidine phosphatase activity. Q9NRX4.

[Learn More](#)

Parameter	Value
Product Name	PHPT1 Rabbit pAb
Remarks/Alias	PHP; PHP14; CGI-202; HSPC141; HEL-S-132P; PHPT1
Species	Human
Gene ID (Human)	29085
Gene ID	29085
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human PHPT1 (NP_054891.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRX4
Protein Weight	14kDa
Shipping	Ice bag