


Anti-Phpt1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nrx4

Item Number: CM0015869



Introduction

Rabbit polyclonal antibody against human PHPT1 (phosphohistidine phosphatase 1). Validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Immunogen: recombinant protein. Research use only.

[Learn More](#)

Parameter	Value
Product Name	PHPT1 Rabbit pAb
Remarks/Alias	PHP; PHP14; CGI-202; HSPC141; HEL-S-132P; PHPT1
Species	Human
Gene ID (Human)	29085
Gene ID	29085
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRX4
Protein Weight	14kDa
Shipping	Ice bag