

Anti-Phrf1 Rabbit Polyclonal Antibody For Wb & Elisa - Q9P1Y6

Item Number: CM0019626

Introduction

Rabbit polyclonal antibody targeting human PHRF1 (PHD and RING finger domain-containing protein 1), validated for Western blot and ELISA. Detects a ~179 kDa band, ideal for studying this metal-binding phosphoprotein.

[Learn More](#)

Parameter	Specification
Product Name	PHRF1 Rabbit pAb
Remarks/Alias	RNF221; PPP1R125; PHRF1
Species	Human
Gene ID (Human)	57661
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1350-1450 of human PHRF1 (NP_001273512.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9P1Y6
Protein Weight	179kDa
Shipping	Ice bag