

Anti-PI3 Kinase P110 Alpha Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P42336

Item Number: CM0007047



Introduction

High-quality rabbit monoclonal antibody against human PI3 Kinase p110 alpha (PIK3CA), validated for WB, IP, ELISA, with cross-reactivity to mouse and rat. Ideal for IVD assay development targeting the PI3K/AKT signaling pathway.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p110 alpha Rabbit mAb
Remarks/Alias	MCM; CCM4; CWS5; MCAP; PI3K; CLAPO; CLOVE; MCMTC; PI3K-alpha; p110-alpha; PI3 Kinase p110 alpha
Species	Human
GeneID (human)	5290
GeneID	5290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 294-610 of human PI3 Kinase p110 alpha (NP_006209.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42336
Protein Weight	124kDa
Shipping	Ice bag