

Anti-Pi3 Kinase P110 Beta Polyclonal Antibody For Wb And Elisa - P42338

Item Number: CM0017387

Introduction

Rabbit polyclonal antibody against human PI3 Kinase p110 beta, suitable for Western blot and ELISA. Targets PIK3CB, involved in cell growth, survival, and oncogenic signaling. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p110 beta Rabbit pAb
Remarks/Alias	PI3K; PIK3C1; P110BETA; PI3KBETA; PI3 Kinase p110 beta
Species	Human
Gene ID (Human)	5291
Gene ID	5291
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-600 of human PI3 Kinase p110 beta (NP_006210.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42338
Protein Weight	123kDa
Shipping	Ice bag