

Anti-PI3 Kinase P110 Beta Rabbit Monoclonal Antibody For Wb And Elisa - P42338

Item Number: CM0006525

Introduction

Anti-PI3 Kinase p110 beta Rabbit Monoclonal Antibody validated for WB and ELISA, detects human, mouse, rat. Ideal for studying PI3K/AKT signaling, cell growth, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p110 beta Rabbit mAb
Remarks/Alias	PI3K; PIK3C1; P110BETA; PI3KBETA; PI3 Kinase p110 beta
Species	Human
GeneID (Human)	5291
GeneID	5291
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-600 of human PI3 Kinase p110 beta (NP_006210.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42338
Protein Weight	123 kDa
Shipping	Ice bag