

Anti-Pi3 Kinase P110 Beta Rabbit Monoclonal Antibody For Wb And Elisa - P42338

Item Number: CM0003802

Introduction

Rabbit monoclonal antibody targeting human PI3 Kinase p110 beta (PIK3CB, P42338) validated for WB and ELISA. Ideal for studying PI3K/AKT signaling, platelet biology, autophagy, and cancer.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p110 beta Rabbit mAb
Remarks/Alias	PI3K; PIK3C1; P110BETA; PI3KBETA; PI3 Kinase p110 beta
Species	Human
GeneID (Human)	5291
GeneID	5291
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PI3 Kinase p110 beta (NP_006210.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42338
Protein Weight	123kDa
Shipping	Ice bag