

Anti-PI3 Kinase P110 Delta Polyclonal Antibody For Wb, Elisa - 000329

Item Number: CM0024282



Introduction

Rabbit polyclonal antibody targeting human PI3 Kinase p110 delta (PIK3CD, 000329). Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for immune signaling research.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p110 delta Rabbit pAb
Remarks/Alias	APDS; PI3K; IMD14; p110D; IMD14A; IMD14B; ROCHIS; P110DELTA; PI3 Kinase p110 delta
Species	Human
GeneID (Human)	5293
GeneID	5293
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-300 of human PI3 Kinase p110 delta (NP_005017.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	000329
Protein Weight	119kDa
Shipping	Ice bag