

Anti-PI3 Kinase P110 Delta Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O00329

Item Number: CM0002979



Introduction

Rabbit monoclonal antibody against human PI3 Kinase p110 delta (O00329), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p110 delta Rabbit mAb
Remarks/Alias	APDS; PI3K; IMD14; p110D; IMD14A; IMD14B; ROCHIS; P110DELTA; PI3 Kinase p110 delta
Species	Human
Gene ID (Human)	5293
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PI3 Kinase p110 delta (O00329).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00329
Protein Weight	119 kDa
Shipping	Ice bag