

Anti-PI3 Kinase P85 Alpha Polyclonal Antibody For Wb, Ihc-P, Elisa - P27986

Item Number: CM0014208



Introduction

Rabbit polyclonal antibody against human PI3 Kinase p85 alpha, validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying PI3K/AKT signaling and insulin metabolic pathways.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 alpha Rabbit pAb
Remarks/Alias	p85; AGM7; GRB1; IMD36; p85alpha; p85-ALPHA; PI3 Kinase p85 alpha
Species	Human
GeneID (Human)	5295
GeneID	5295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-380 of human PI3 Kinase p85 alpha (NP_852556.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27986
Protein Weight	84kDa
Shipping	Ice bag