

Anti-Pi3 Kinase P85 Alpha Rabbit Mab For Wb, If, Ihc, Elisa - P27986

Item Number: CM0005565



Introduction

High-quality recombinant PI3 Kinase p85 alpha rabbit monoclonal antibody, validated for WB, IF, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Suitable for PI3K/AKT signaling and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 alpha Rabbit mAb
Remarks/Alias	p85; AGM7; GRB1; IMD36; p85alpha; p85-ALPHA; PI3 Kinase p85 alpha
Species	Human
Gene ID	5295 (Human)
Immunogen	Recombinant fusion protein containing amino acids 240-380 of human PI3 Kinase p85 alpha (NP_852664.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27986
Protein Weight	84kDa
Shipping	Ice bag