

Anti-PI3 Kinase P85 Alpha Rabbit Monoclonal Antibody For Wb And Elisa - P27986

Item Number: CM0005199



Introduction

Monoclonal rabbit antibody targeting PI3 Kinase p85 alpha, validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying PI3K/AKT signaling, insulin response, and receptor tyrosine kinase pathways.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 alpha Rabbit mAb
Remarks/Alias	p85; AGM7; GRB1; IMD36; p85alpha; p85-ALPHA
Species	Human
Gene ID (Human)	5295
Gene ID	5295
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27986
Protein Weight	84kDa
Shipping	Ice bag