

Anti-Pi3 Kinase P85 Alpha/Beta Rabbit Monoclonal Antibody For Wb, Elisa - P27986 / O00459

Item Number: CM0007613



Introduction

High-quality rabbit monoclonal antibody targeting PI3 Kinase p85 alpha and beta (PIK3R1/PIK3R2) for Western blot and ELISA applications. Validated in human and mouse species. Ideal for PI3K/AKT pathway research.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 alpha/beta Rabbit mAb
Remarks/Alias	PIK3R1/PIK3R2; PI3 Kinase p85 alpha/beta
Species	Human
Gene ID (Human)	5295/5296
Gene ID	5295/5296
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P27986/O00459
Protein Weight	42-54 kDa/81-84 kDa
Shipping	Ice bag