


Anti-Pi3 Kinase P85 Beta Monoclonal Antibody For Wb, Elisa - O00459

Item Number: CM0012753



Introduction

Rabbit monoclonal antibody against PI3 Kinase p85 beta, validated in WB and ELISA. Detects human and rat p85 β (82 kDa). Ideal for research on PI3K/AKT signaling, cancer, and metabolic pathways.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 beta Rabbit mAb
Remarks/Alias	p85; MPPH; P85B; MPPH1; p85beta; p85-BETA; PI3 Kinase p85 beta
Species	Human
GeneID (Human)	5296
GeneID	5296
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human PI3 Kinase p85 beta (O00459).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O00459
Protein Weight	82kDa
Shipping	Ice bag