


Anti-PI3 Kinase P85 Beta Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O00459

Item Number: CM0021029



Introduction

Rabbit polyclonal antibody against PI3 Kinase p85 beta (PIK3R2), validated for WB, IHC-P, ELISA in Human, Mouse, Rat. Regulatory subunit of PI3K, involved in cell growth and signaling. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 beta Rabbit pAb
Remarks/Alias	p85; MPPH; P85B; MPPH1; p85beta; p85-BETA; PI3 Kinase p85 beta
Species	Human
GeneID (Human)	5296
GeneID	5296
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 108-298 of human PI3 Kinase p85 beta (NP_005018.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00459
Protein Weight	82kDa
Shipping	Ice bag