


Anti-Pi3 Kinase P85 Beta Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00459

Item Number: CM0021714



Introduction

Rabbit polyclonal antibody targeting human PI3 Kinase p85 beta, validated for WB, IHC-P, IF/ICC, ELISA. Detects endogenous p85 beta in human, mouse, and rat samples. Suitable for PI3K/AKT signaling and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 beta Rabbit pAb
Remarks/Alias	p85; MPPH; P85B; MPPH1; p85beta; p85-BETA; PI3 Kinase p85 beta
Species	Human
Gene ID (Human)	5296
Gene ID	5296
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human PI3 Kinase p85 beta (NP_005018.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00459
Protein Weight	82kDa
Shipping	Ice bag