


Anti-PI3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P19957

Item Number: CM0028262



Introduction

Rabbit polyclonal antibody against human PI3 (Elafin, SKALP, ESI). Validated for Western blot, immunohistochemistry, and ELISA. Cross-reacts with mouse. Ideal for investigating elastase inhibition and inflammatory skin diseases.

[Learn More](#)

Parameter	Value
Product Name	PI3 Rabbit pAb
Remarks/Alias	ESI; WAP3; SKALP; WFDC14; cementoin; PI3
Species	Human
Gene ID	5266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-117 of human PI3 (NP_002629.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19957
Protein Weight	12kDa
Shipping	Ice bag