

Anti-Pi3K-Gamma Rabbit Monoclonal Antibody For Wb And Elisa - P48736

Item Number: CM0004262



Introduction

High-quality rabbit monoclonal antibody against PI3K-gamma (p110 γ , 126kDa) for WB and ELISA applications. Detects human and mouse PI3K γ ; ideal for studying immune cell signaling, chemotaxis, and GPCR-mediated pathways.

[Learn More](#)

Parameter	Specification
Product Name	PI3K-gamma Rabbit mAb
Remarks/Alias	PI3K; PIK3; IMD97; PI3CG; PI3Kgamma; p110gamma; p120-PI3K
Species	Human
GeneID (Human)	5294
GeneID	5294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-342 of human PI3 Kinase-gamma (NP_002640.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48736
Protein Weight	126kDa
Shipping	Ice bag