

Anti-Pi3K P85A/P55γ/P85B Polyclonal Antibody For Wb, If/Icc, Elisa - P27986 / O00459 / Q92569

Item Number: CM0023320



Introduction

Rabbit polyclonal antibody targeting PI3K regulatory subunits p85α, p55γ, and p85β (PIK3R1, PIK3R2, PIK3R3). Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, and rat. Ideal for insulin signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PI3KP85α/P55γ/P85β Rabbit pAb
Remarks/Alias	PIK3R1/PIK3R2/PIK3R3; PI3KP85α/P55γ/P85β
Species	Human
Gene ID (Human)	5295/5296/8503
Gene ID	8503/5295/5296
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400 to 500 of human PIK3R1/PIK3R2/PIK3R3 (NP_852664.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27986/O00459/Q92569
Protein Weight	42-54 kDa/81-84 kDa
Shipping	Ice bag