


Anti-PI4K2B Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tcg2

Item Number: CM0023528



Introduction

Rabbit polyclonal antibody targeting human PI4K2B, suitable for Western blot and ELISA. Recognizes phosphatidylinositol 4-kinase type 2-beta, involved in PI4P synthesis and vesicular trafficking.

[Learn More](#)

Parameter	Value
Product Name	PI4K2B Rabbit pAb
Remarks/Alias	PIK42B; PI4KIIB; PI4K2B
Species	Human
Gene ID (Human)	55300
Gene ID	55300
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human PI4K2B (NP_060793.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TCG2
Protein Weight	55kDa
Shipping	Ice bag