


Anti-Pi4Ka Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P42356

Item Number: CM0021971



Introduction

Rabbit polyclonal antibody targeting human PI4KA, validated for Western blot, immunofluorescence/ICC, and ELISA. Detects human, mouse, and rat PI4KA. Ideal for phosphatidylinositol signaling and spastic paraplegia research.

[Learn More](#)

Parameter	Value
Product Name	PI4KA Rabbit pAb
Remarks/Alias	SPG84; GIDID2; PIK4CA; PMGYCHA; pi4K230; PI4K-ALPHA; PI4KA
Species	Human
Gene ID (Human)	5297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1780-1860 of human PI4KA (P42356)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42356
Protein Weight	237kDa
Shipping	Ice bag